

# Esy-Pik Microtube Picker

The Vitl Esy-Pik  
Microtube Picker allows  
one-handed removal of  
microtubes from a rack




## Description

The Vitl Esy-Pik microtube picker allows one-handed removal of microtubes and incorporates an eject button to allow pick and place operations whilst avoiding contact with the tube.

This minimizes the risk of heating a frozen sample and possible cross contamination.

## Benefits

- Single-handed operation
- Temperature and chemically resistant tube contact area
- Tube selection from anywhere in the rack
- Reduced risk of thawing of frozen tubes and sample contamination
- Easy-grip
- Lightweight
- Can be used with laboratory gloves



## Technical Specifications

### Weight

**MicroTube Esy-Pik** 20g

**Cryotube Esy-Pik** 23g

### Dimensions (mm) L x Dia. Excluding pocket clip

**MicroTube Esy-Pik** 150 x 17

**Cryotube Esy-Pik** 165 x 17

## Catalogue number

**V201001** Vitl Microtube Esy-Pik (for 8.5mm tubes)

**V201002** Vitl Cryotube Esy-Pik (for Cryotubes up to 12.5mm cap)

Please contact [sales@vitlproducts.com](mailto:sales@vitlproducts.com)  
for ordering details for the Esy-Pik

**Vitl**   
LIFE SCIENCE SOLUTIONS  
A member of the ITL group

T: +44 1233 638 383 **UK**  
T: +1 804 381 0905 **USA**

[vitlproducts.com](http://vitlproducts.com)