

	Peelable?	Pierceable?	Temperature range	Seals Polypropylene?	Seals Polyethylene?	Seals Polystyrene?	Seals COC?	Color code	Catalogue Number	Product Notes
20µm Heat Sealing Foil	X	✓	-20 - +120°C	✓	✓	✓	X	Green	V901002	
38µm Heat Sealing Foil	✓	✓	-20 - +110°C	✓	✓	✓	X	Blue	V901003	Only peels from PS
70µm Heat Sealing Foil	✓	X	-80 - +80°C	✓	✓	X	X	Red	V901001	DMSO resistant
75µm Clear Sealing Film	X	✓	-80 - +80°C	✓	✓	✓	✓	75µm	V901005	Suitable for PCR at 110°C
85µm Heat Sealing Foil	✓	X	-200 - +110°C	✓	X	X	X	Black	V901004	
105µm Clear Sealing Film	✓	X	-80 - +180°C	✓	✓	✓	✓	105µm	V901006	Suitable for PCR at 110°C