

Heat Sealing Foils & Films

Product	Peelable	Piercable	Temp Range	Polyethylene Plate (PE)	Polypropylene Plate (PP)	Polystyrene Plate (PS)	COC	Colour Code	Catalogue Number
20µm Heat Sealing Foil	X	✓	-20 to +120°C	✓	✓	✓	X	Green	V901002
38µm Heat Sealing Foil	✓	✓	-20 to +110°C	✓	✓	✓	X	Blue	V901003
68µm Heat Sealing Foil	✓	X	-80 to +110°C	✓	✓	✓	✓	Yellow	V901008
70µm Heat Sealing Foil	✓	✓	-80 to +80°C	✓	✓	✓	✓	Orange	V901010
70µm Heat Sealing Foil	✓	X	-80 to +80°C	✓	✓	X	X	Red	V901001
75µm Clear Sealing Film	X	✓	-80 to +80°C	✓	✓	✓	✓	75µm	V901005
85µm Heat Sealing Foil	✓	X	-200 to +110°C	✓	X	X	X	Black	V901004
105µm Clear Sealing Film	✓	X	-80 to +180°C	✓	✓	✓	✓	105µm	V901006

*All dimensions are nominal. Some variation to be expected.

V901020 Foil Weight (sold separately)

Not sure which heat seal is right for you? Contact us for a free sample pack

Tel: +44 (0)1233 638383
 Email: sales@vitlproducts.com
 Website: www.vitlproducts.com

